

Investigations on human biomonitoring for assessing the exposure to synthetic rose oxide in the general population

C. Tiwald, V. Spindler, M. Scherer, N. Pluym, O. Peschel, E. Leibold, G. Scherer

Analytisch-Biologisches Forschungslabor GmbH, Semmelweisstr. 5, 82152 Planegg, Germany

ISBM 13, Milan 2025, 10-September-2025

Rose Oxide (RO)

Use as a fragrance in cosmetic, cleaning agents, household products, etc.

NATURAL VS. SYNTHETIC

mainly (-)-cis-isomer

↓
0.3 – 1.2 % natural rose oxide in roses → time-consuming and expensive production


equal parts
(-)- and (+)-cis-isomers
& (-)- and (+)-trans isomers

↓
(+)-isomers specific for synthetic rose oxide

TOXICITY

- Low acute toxicity ($LD_{50} > 2000$ mg/kg body weight)
- ECHA: NOAEL 100 mg/kg bw/d (♂), 300 mg/kg bw/d (♀) (oral subacute dosage)
- Skin and eye irritating, not skin-sensitizing, no mutagenic effects, may impair fertility

MW: 154 g/mol



2-(2-Methyl-1-propenyl)-4-methyl-tetrahydropyran

High production chemical: Broad exposure of the general population expected

➔ **HUMAN BIOMONITORING**


Workflow

Identification of a suitable metabolite

Method development and validation

Excretion kinetics after single oral and dermal application

Application to urine samples



Rose oxide

Postulated Metabolism


- Postulated from in vitro experiments (microsomes from mouse, rat, rabbit & human)
- 10 postulated metabolites
 - Hydroxylation and oxidation reactions
 - Incubations with human microsomes: **2, 3, 4, 5, 11**
 ➔ **2, 3, 5** and **11** retain both stereocenters

Identification of a suitable metabolite

Method development and validation

Excretion kinetics

Application to urine samples



Postulated Metabolism


- Postulated from in vitro experiments (microsomes from mouse, rat, rabbit & human)
- 10 postulated metabolites
 - Hydroxylation and oxidation reactions
 - Incubations with human microsomes:
2, 3, 4, 5, 11
 ➡ **2, 3, 5** and **11** retain both stereocenters

Identification of a suitable metabolite

Method development and validation

Excretion kinetics

Application to urine samples



In-house Metabolism Study

- 5 healthy volunteers dosed with (+)-RO (oral and dermal)
- Urine collection for 48 h (oral) and 72 h (dermal)
Documentation of time points and urine volume
- Non-targeted analysis of selected urine fractions with UPLC-HRMS
- Data evaluation – **Screening for metabolites**

Identification of a suitable metabolite

Method development and validation

Excretion kinetics

Application to urine samples

oral
0.1 mg/kg bw



dermal
0.5 mg/kg bw



Non-targeted analysis



Q-Exactive-HF-X Orbitrap MS

Screening for suitable metabolites

Full MS1 and MS2 spectrum of selected candidates

Results – Non-targeted Analysis




- 18 features followed the typical excretion pattern – 4 metabolites
- m/z 201: “Hydroxy-Carboxy-RO”, [C₁₀H₁₆O₄], [**5**]
- m/z 153: “9-OH-RO” –H₂O, [C₁₀H₁₆O], [**11**]
- m/z 169: “Di-OH-RO”/Lactone, [C₁₀H₁₆O₂], [**9/2''/2'**]
- m/z 229: not identified

Identification of a suitable metabolite

Method development and validation

Excretion kinetics

Application to urine samples



Metabolite Selection

Criteria for suitable metabolites

Clarified structure
Confirmed by reference standard



Both stereocenters remain intact
((+)-cis-isomer is specific for synthetic rose oxide)




Metabolite formation sufficient for analysis

Identification of a suitable
metabolite

Method
development and
validation

Excretion kinetics


Application to
urine samples



×

✓


?



×

✓


?



×

×

?



✓

✓

×



Parent compound (+)-cis-RO

Quantitative GC-MS Method

Analytical Workflow

Sample workup


3 mL urine sample
+ 1 g KCl
+ IS (linalool-d₅)
(10 mL Headspace-Vial)



SPME-HS-GC-MS

EI, SIM

Rose oxide (Quan)	m/z: 139
Rose oxide (Qual)	m/z: 69
Linalool-d ₅	m/z: 159



Identification of a
suitable metabolite

Method development and
validation

Excretion kinetics

Application to
urine samples


Method Validation

According to FDA Guideline on Bioanalytical Method Validation (2001)

Calibration range	1 – 200 ng/L		
LLOQ	1 ng/L		
Accuracy and precision	level	accuracy [%]	precision (CV, [%])
Inter-day (N = 3x5)	1 ng/L	99.7	11.2
	2 ng/L	91.2	8.7
	10 ng/L	96.2	5.6
	100 ng/L	99.9	4.2

Parameters: Linearity, accuracy & precision, selectivity, matrix effect, stability, carryover, dilution integrity

Urine sample,
200 ng/L



In-house Metabolism Study

- 5 healthy volunteers dosed with (+)-RO (oral and dermal)
- Urine collection for 48 h (oral) and 72 h (dermal)
Documentation of time points and urine volume
- Targeted analysis of all urine fractions with SPME-HS-GC-MS
- Data evaluation – **Excretion kinetics**

oral
0.1 mg/kg bw



dermal
0.5 mg/kg bw



Targeted
analysis






Agilent GCMS 5975

Evaluation of excretion kinetics and
toxicokinetic parameters for each
participant

Excretion Kinetics - Oral Application

- Baseline levels of (+)-cis-RO before application: 0 – 6.2 ng/L
- Fast elimination
 - Mean $t_{\max}$ = 1.6 h
 - Complete elimination after 48 h
- Recovery of 7 ppm
- Bi-exponential decline with 2 elimination phases



	Mean $\pm$ SD	Median	Min–Max
Cumulative excretion after 48 h [pmol] (A_{48h})	337 $\pm$ 196	296	71–657
$t_{\max}$ [h]	1.63 $\pm$ 0.98	1.25	0.73–3.33
Total (+)-cis-RO excreted after 3 h [%]	62 $\pm$ 23	67	22–89
Total (+)-cis-RO excreted after 6 h [%]	76 $\pm$ 17	86	44–92
Total (+)-cis-RO excreted after 12 h [%]	83 $\pm$ 1	91	55–95
Total (+)-cis-RO excreted after 24 h [%]	90 $\pm$ 8	94	77–98
Total (+)-cis-RO after 48 h [%]	99 $\pm$ 2	100	96–100
Elimination constant λ_1 [h ⁻¹]	0.36 $\pm$ 0.19	0.35	0.10–0.69
Elimination half-life $t_{1/2 \lambda_1}$ [h]	4.10 $\pm$ 3.01	2.88	1.44–9.93
Elimination constant λ_2 [h ⁻¹]	0.07 $\pm$ 0.02	0.06	0.05–0.09
Elimination half-life $t_{1/2 \lambda_2}$ [h]	14.7 $\pm$ 3.35	15.9	10.6–18.0
Urinary excretion factor F_{ue} (24 h) [%]	0.007 $\pm$ 0.005	0.005	0.001–0.017
Urinary excretion factor F_{ue} (48 h) [%]	0.007 $\pm$ 0.005	0.005	0.001–0.017

Excretion Kinetics - Dermal Application

- Baseline levels of (+)-cis-RO application: 0.7 – 4.8 ng/L
- Slower elimination compared to oral application
 - Mean $t_{\max}$ = 2.77 h
 - 82 % of (+)-cis-RO eliminated after 48 h
- Recovery of 0.3 ppm
- Dermal resorption rate: 18 %



Participant	Dermal (+)-RO dose [μmol]	Excreted amount after dermal application [nmol]	Excreted amount after oral application [%]	Resorption rate [%]
1	259.3	0.77	0.0008	37.1
2	259.3	0.19	0.0005	14.7
3	275.5	0.10	0.0005	7.3
4	201.0	0.08	0.0017	2.3
5	243.1	0.07	0.0001	28.8
Mean $\pm$ SD				18.0% $\pm$ 13.1%

Identification of a suitable metabolite

Method development and validation

Excretion kinetics

Application to urine samples



	Mean $\pm$ SD	Median	Min–Max
Cumulative excretion after 72 h [pmol] (A_{72h})	103 $\pm$ 44.6	83.0	65–189
$t_{\max}$ [h]	2.77 $\pm$ 0.73	3.00	1.67–3.50
Total (+)-cis-RO excreted after 3 h [%]	29 $\pm$ 12	29	13–45
Total (+)-cis-RO excreted after 6 h [%]	35 $\pm$ 10	30	22–50
Total (+)-cis-RO excreted after 12 h [%]	50 $\pm$ 12	47	37–73
Total (+)-cis-RO excreted after 24 h [%]	63 $\pm$ 8	60	57–78
Total (+)-cis-RO excreted after 48 h [%]	82 $\pm$ 7	86	72–90
Urinary excretion factor F_{ue} (24 h) [‰]	0.0003 $\pm$ 0.0001	0.0002	0.0002–0.0004
Urinary excretion factor F_{ue} (48 h) [‰]	0.0003 $\pm$ 0.0001	0.0003	0.0002–0.0006

Application

- “General population”

- 13 non-occupationally exposed participants
- 1 participant intentionally exposed to rose oxide-containing perfume (positive control)

- Quantifiable amounts of (+)-cis-RO in 54 % of unexposed samples

- 1.2 ng/L on average (LLOQ: 1 ng/L)



positive control




SPME-HS-GC-MS
Targeted analysis



**Very low urinary excretion rate of
7 ppm (oral) / 0.3 ppm (dermal)**

Summary and Conclusion



CHALLENGES FOR HBM

Low excretion rate of (+)-cis-RO
Possible contamination by environmental RO

Acknowledgment

- Funding

This study was funded by the Chemie Wirtschaftsförderungsgesellschaft, Frankfurt/Main, Germany. The method was developed as part of a 10-year human biomonitoring initiative between the Federal Ministry for Environment, Nature Conservation and Nuclear Safety (BMUV) and the Verband der chemischen Industrie e.V. (German Chemical Industry Association (VCI)).

- ABF Team

Veronika Spindler
Dr. Therese Burkhardt
Dr. Nikola Pluym
Prof. Dr. Gerhard Scherer
Dr. Max Scherer

